

# BASELINE Acetic Acid MS      Maximum Specifications and Typical Values

<b>SEASTAR™ BASELINE®</b>	
Product Code	S020601
Assay (CH <sub>3</sub> COOH, w/w):	≥ 99%
Analyte	Trace Impurities in ppt (pg/g)
Aluminum (Al)	50
Antimony (Sb)	50
Arsenic (As)	50
Barium (Ba)	10
Beryllium (Be)	10
Bismuth (Bi)	10
Cadmium (Cd)	10
Calcium (Ca)	50
Cerium (Ce)	10
Cesium (Cs)	10
Chromium (Cr)	10
Cobalt (Co)	10
Copper (Cu)	10
Dysprosium (Dy)	1
Erbium (Er)	1
Europium (Eu)	1
Gadolinium (Gd)	1
Gallium (Ga)	10
Germanium (Ge)	10
Hafnium (Hf)	10
Holmium (Ho)	1
Indium (In)	1
Iron (Fe)	50
Lanthanum (La)	1
Lead (Pb)	10
Lithium (Li)	10
Lutetium (Lu)	10
Magnesium (Mg)	50
Manganese (Mn)	10
Molybdenum (Mo)	10
Neodymium (Nd)	1
Nickel (Ni)	50
Platinum (Pt)	50
Potassium (K)	50
Praseodymium (Pr)	1
Rhenium (Re)	10
Rhodium (Rh)	50
Rubidium (Rb)	10
Ruthenium (Ru)	50
Samarium (Sm)	1
Scandium (Sc)	10
Selenium (Se)	Information only
Silver (Ag)	50
Sodium (Na)	100
Strontium (Sr)	10
Tellurium (Te)	1
Terbium (Tb)	1
Thallium (Tl)	10
Thorium (Th)	1
Thulium (Tm)	1
Tin (Sn)	50
Titanium (Ti)	10
Tungsten (W)	10
Uranium (U)	1
Vanadium (V)	10
Ytterbium (Yb)	1
Yttrium (Y)	1
Zinc (Zn)	50
Zirconium (Zr)	10