

BASELINE Hydrogen Peroxide Maximum Specifications

SEASTAR™ BASELINE®	
Product Code	S021001
Assay (H ₂ O ₂ w/w):	30 - 32%
Analyte	Trace Impurities in ppt (pg/g)
Aluminum (Al)	50
Antimony (Sb)	10
Arsenic (As)	100
Barium (Ba)	10
Beryllium (Be)	10
Bismuth (Bi)	10
Boron (B)	100
Cadmium (Cd)	10
Calcium (Ca)	100
Cerium (Ce)	1
Cesium (Cs)	1
Chromium (Cr)	10
Cobalt (Co)	10
Copper (Cu)	10
Dysprosium (Dy)	1
Erbium (Er)	1
Europium (Eu)	1
Gadolinium (Gd)	1
Gallium (Ga)	10
Germanium (Ge)	10
Gold (Au)	10
Hafnium (Hf)	1
Holmium (Ho)	1
Indium (In)	1
Iron (Fe)	20
Lanthanum (La)	1
Lead (Pb)	10
Lithium (Li)	10
Lutetium (Lu)	1
Magnesium (Mg)	20
Manganese (Mn)	10
Mercury (Hg)	50
Molybdenum (Mo)	10
Neodymium (Nd)	1
Nickel (Ni)	20
Niobium (Nb)	10
Palladium (Pd)	10
Potassium (K)	20
Praseodymium (Pr)	1
Rhenium (Re)	10
Rhodium (Rh)	10
Rubidium (Rb)	10
Ruthenium (Ru)	10
Samarium (Sm)	1
Scandium (Sc)	10
Selenium (Se)	100
Silver (Ag)	10
Sodium (Na)	50
Strontium (Sr)	10
Tantalum (Ta)	10
Tellurium (Te)	1
Terbium (Tb)	1
Thallium (Tl)	1
Thorium (Th)	1
Thulium (Tm)	1
Tin (Sn)	50
Titanium (Ti)	20
Tungsten (W)	20
Uranium (U)	1
Vanadium (V)	10
Ytterbium (Yb)	1
Yttrium (Y)	1
Zinc (Zn)	50
Zirconium (Zr)	10