

BASELINE Sulphuric Acid MS Maximum Specifications

SEASTAR™ BASELINE®	
Product Code	S020301
Assay (H ₂ SO ₄ w/w):	93 - 98%
Analyte	Trace Impurities in ppt (pg/g)
Aluminum (Al)	50
Antimony (Sb)	50
Arsenic (As)	500
Barium (Ba)	10
Beryllium (Be)	10
Bismuth (Bi)	10
Cadmium (Cd)	10
Calcium (Ca)	50
Cerium (Ce)	10
Cesium (Cs)	10
Chromium (Cr)	10
Cobalt (Co)	10
Copper (Cu)	10
Dysprosium (Dy)	10
Erbium (Er)	10
Europium (Eu)	10
Gadolinium (Gd)	10
Gallium (Ga)	10
Germanium (Ge)	100
Hafnium (Hf)	10
Holmium (Ho)	10
Indium (In)	10
Iron (Fe)	50
Lanthanum (La)	10
Lead (Pb)	10
Lithium (Li)	10
Lutetium (Lu)	10
Magnesium (Mg)	50
Manganese (Mn)	10
Mercury (Hg)	100
Molybdenum (Mo)	10
Neodymium (Nd)	10
Nickel (Ni)	50
Niobium (Nb)	10
Palladium (Pd)	Information Only
Platinum (Pt)	Information Only
Potassium (K)	50
Praseodymium (Pr)	10
Rhodium (Rh)	50
Rubidium (Rb)	10
Samarium (Sm)	10
Scandium (Sc)	10
Selenium (Se)	500
Silver (Ag)	50
Sodium (Na)	50
Strontium (Sr)	10
Tantalum (Ta)	Information Only
Tellurium (Te)	50
Terbium (Tb)	10
Thallium (Tl)	10
Thorium (Th)	10
Thulium (Tm)	10
Tin (Sn)	50
Titanium (Ti)	50
Tungsten (W)	10
Uranium (U)	10
Vanadium (V)	10
Ytterbium (Yb)	10
Yttrium (Y)	10
Zinc (Zn)	50
Zirconium (Zr)	10